DCM074-5
Ed. 01/2015

S100B ELISA

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11. Rosén H et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicenter study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Stigbrand T, Nyberg L, Ullen A, Haglid K, Sandstrom E, Brundell J. A new specific method for measuring S-100B in serum. *Int J Biol Markers*. 2000 Jan-Mar;15(1):33-40.
15. Chen DQ, Zhu LL. Dynamic change of serum protein S100b and its clinical significance in patients with traumatic brain injury.. *Chin J Traumatol*, 2005 Aug 1;8(4):245-8.
16. Sawauchi S, Taya K, Murakami S, Ishi T, Ohtsuka T, Kato N, Kaku S, Tanaka T, Morooka S, Yuhki K, Urashima M, Abe T. Serum S-100B protein and neuron-specific enolase after traumatic brain injury. *No Shinkei Geka*. 2005 Nov;33(11):1073-80

BIBLIOGRAFIA

1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins: structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E, et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.

Ed. 01/2015

DCM074-5

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM074-5
Ed. 01/2015

S100B ELISA

for routine analysis

Direct immunoenzymatic determination of S100B in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO074

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of S100B concentration in human serum or plasma.

S100B kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

S100 is a 20 kDa protein belonging to the S100/calmodulin/troponin C superfamily of EF-hand calcium-binding proteins. S100 was originally isolated from human brain and considered a glial-cell specific protein (1). Today, 20 monomers of the S100 family have been identified based on structural and functional similarities (2,3). Most of the S100 proteins exist as dimers and are expressed in a cell-specific manner. Two of the S100 monomers, designated S100α and S100β (4) are highly conserved between species and are found as homo- (ββ) and heterodimers (αβ) in central nervous system glial cells and in certain peripheral cells eg. Schwann cells, melanocytes, adipocytes, and chondrocytes (5). S100αβ and S100ββ are also present in malignant tissues, most notably in melanoma and to a lesser extent in glioma, thyroid cell carcinoma and renal cell carcinoma (2). Determination of S100B (like S100αβ and S100ββ units) in serum has been shown to be clinically useful for prognosis and treatment monitoring of patients diagnosed with malignant melanoma (6-9). Studies also suggest that S100B may be useful in the management of patients with brain damage from eg. traumatic head injury, perinatal asphyxia, cardiac arrest, cardiac surgery and stroke (10-13).

2. PRINCIPLE

S100B ELISA test is based on binding of S100B by two antibodies, one immobilized on microwell plates, and the other one conjugated with horseradish peroxidase (HRP). The assay is on two steps binding procedure and, after every incubation step, the bound/free separation is performed by a simple solid-phase washing. Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution (H₂SO₄) is added. The color intensity is proportional to the S100B concentration in the sample. S100B concentration in the sample is calculated through a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1 Reagents and materials supplied in the kit.

- Calibrators** (6 vials, lyophilized)

CAL0	REF DCE002/7406-0
CAL1	REF DCE002/7407-0
CAL2	REF DCE002/7408-0
CAL3	REF DCE002/7409-0
CAL4	REF DCE002/7410-0
CAL5	REF DCE002/7411-0
- Controls** (2 vials, lyophilized)

Negative Control	REF DCE045/7401-0
Positive Control	REF DCE045/7402-0
- Conjugate buffer** (1 vial, 20 mL)
Tris buffer; BSA 10 g/L, tween 0.05%
REF DCE044-0
- Conjugate** (1 vial, 1 mL)
Anti S100B antibody conjugated with horseradish peroxidase (HRP)
REF DCE002/7402-0
- Assay buffer** (1 vial, 12 mL)
Tris buffer; BSA 10 g/L
REF DCE043-0
- Coated Microplate** (1 breakable microplate)
Anti S100B antibody adsorbed on microplate
REF DCE002/7403-0
- TMB Substrate** (1 vial, 15 mL)
H₂O₂-TMB 0.26 g/L (avoid any skin contact)
REF DCE004-0
- Stop Solution** (1 vial, 15 mL)
Sulphuric acid 0.15 mol/L (avoid any skin contact)
REF DCE005-0
- 50X Conc. Wash Solution** (1 vial, 20 mL)
NaCl 45g/L Tween-20 55g/L
REF DCE006-0

3.2 Necessary reagents not supplied with the kit

Distilled water.

3.3 Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Notes

Store all reagents between 2÷8°C in the dark.

Open the bag of reagent 6 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened it is stable at 2-8°C until the expiry date of the kit.

Do not remove the adhesive sheets on the strips unutilized.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Sodium Merthiolate or Proclin 300^R as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of S100B from 10 to 5000 pg/mL.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If

more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1 Preparation of the sample

The S100B determination can be carried out in human serum or plasma. Do not use hemolyzed samples.

Samples can be stored at 2÷8°C for 1 day; for long periods store at -20°C. Avoid repetitive freezing and thawing of samples. Do not leave the samples at room temperature (22-28°C) for long period.

For sample with concentration over 5 ng/mL dilute the sample with Assay buffer

6.2 Preparation of the Calibrators and Controls

Reconstitute the Calibrators and Controls with 1 mL of distilled water before use.

Once reconstituted they are stable for 4 weeks at 2÷8°C and about six months if stored at -20°C. It is advised to aliquot the reconstituted content in small aliquots and store them at -20°C. Avoid repetitive freezing and thawing and long time exposure at room temperature (22-28°C).

The value of Calibrators concentration is lot specific and is reported on the Calibrators vial labels and on the Certificate of Analysis.

6.3 Preparation of the Conjugate

Prepare 2 hours before use.

Add 50 µL of Conjugate (reagent 4) to 1 mL of Conjugate Buffer (reagent 3). The quantity of diluted conjugate is proportional at the number of the tests. Mix gently for 5 minutes, with a rotating mixer. Stable for 3 hours at room temperature (22-28°C).

6.4 Preparation of the wash solution

Dilute the content of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C.

6.5 Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₅	50 µL		
Sample/ Controls		50 µL	
Assay Buffer	50 µL	50 µL	
Incubate 2 h at room temperature (22±28°C). Remove the content from each well, wash the wells 6 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
Diluted Conjugate	100 µL	100 µL	
Incubate 1 h at room temperature (22±28°C). Remove the content from each well, wash the wells 6 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 30 minutes at room temperature (22±28°C), in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. RESULTS

7.1 Mean Absorbance

Calculate the mean of the absorbencies (E_m) for each point of the calibration curve (C₀-C₅) and of each sample.

7.2 Calibration curve

Plot the values of absorbance (E_m) of the Calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points (ex: Cubic spline or Four Parameter Logistic).

7.3 Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations in pg/mL.

8. REFERENCES VALUES

Each laboratory must establish its own normal ranges based on patient population:

Normal range: Mean = 50 pg/mL SD = 15 pg/mL

Pathological limit: >75 pg/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

9. PERFORMANCE CHARACTERISTICS

9.1 Specificity

The antibody recognizes specifically the β subunit therefore reacts with S100αβ and S100ββ units, it is not reactive against the S100αα unit. Cross-reacts with S100 from bovine, porcine, rabbit cat and rat. Does not reacts with the other EF-hand family proteins.

9.2 Sensitivity

The lowest detectable concentration of S100B that can be distinguished from the Calibrator 0 is 35,27 pg/mL.

9.3 Hook Effect

In this assay, no hook effect is observed up to 5000 pg/mL of S100B

10. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local

BIBLIOGRAPHY

1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins: structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E, et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicenter study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Stigbrand T, Nyberg L, Ullen A, Haglid K, Sandstrom E, Brundell J. A new specific method for measuring S-100B in serum. *Int J Biol Markers*. 2000 Jan-Mar;15(1):33-40.
15. Chen DQ, Zhu LL. Dynamic change of serum protein S100b and its clinical significance in patients with traumatic brain injury.. *Chin J Traumatol*, 2005 Aug 1;8(4):245-8.
16. Sawauchi S, Taya K, Murakami S, Ishi T, Ohtsuka T, Kato N, Kaku S, Tanaka T, Morooka S, Yuhki K, Urashima M, Abe T. Serum S-100B protein and neuron-specific enolase after traumatic brain injury. *No Shinkei Geka*. 2005 Nov;33(11):1073-80

Ed. 01/2015

DCM074-5

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

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1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins: structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.

6. Banfalvi T et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E, et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicenter study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Stigbrand T, Nyberg L, Ullen A, Haglid K, Sandstrom E, Brundell J. A new specific method for measuring S-100B in serum. *Int J Biol Markers*. 2000 Jan-Mar;15(1):33-40.
15. Chen DQ, Zhu LL. Dynamic change of serum protein S100b and its clinical significance in patients with traumatic brain injury.. *Chin J Traumatol*, 2005 Aug 1;8(4):245-8.
16. Sawauchi S, Taya K, Murakami S, Ishi T, Ohtsuka T, Kato N, Kaku S, Tanaka T, Morooka S, Yuhki K, Urashima M, Abe T. Serum S-100B protein and neuron-specific enolase after traumatic brain injury. *No Shinkei Geka*. 2005 Nov;33(11):1073-80

Ed. 01/2015

DCM074-5

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

	DE ES FR GB IT PT	In vitro Diagnostikum Producto sanitario para diagnóstico In vitro Dispositif medical de diagnostic in vitro In vitro Diagnostic Medical Device Dispositivo medico-diagnostico in vitro Dispositivos medicos de diagnostico in vitro		DE ES FR GB IT PT	Hergestellt von Elaborado por Fabriqué par Manufacturer Produttore Produzido por
	DE ES FR GB IT PT	Achtung, Begleitdokumente Precaución, consulte los documentos adjuntos Attention, veuillez consulter les documents d'accompagnement Caution, consult accompanying documents Attenzione, consultare la documentazione allegata Atenção, consultar os documentos de acompanhamento	 yyyy-mm	DE ES FR GB IT PT	Herstellungs datum Fecha de fabricacion Date de fabrication Date of manufacture Data di produzione Data de produção
 yyyy-mm-dd	DE ES FR GB IT PT	Verwendbar bis Estable hasta (usar antes de último día del mes) Utiliser avant (dernier jour du mois indiqué) Use by (last day of the month) Utilizzare prima del (ultimo giorno del mese) Utilizar (antes ultimo dia do mês)		DE ES FR GB IT PT	Biogefährdung Riesco biológico Risque biologique Biological risk Rischio biologico Risco biológico
	DE ES FR GB IT PT	Gebrauchsanweisung beachten Consultar las instrucciones Consulter le mode d'emploi Consult instructions for use Consultare le istruzioni per l'uso Consultar instruções para uso		DE ES FR GB IT PT	Chargenbezeichnung Codigo de lote Numero de lot Batch code Codice del lotto Codigo do lote
 $\Sigma = xx$	DE ES FR GB IT PT	Ausreichend für "n" Tests Contenido suficiente para "n" tests Contenu suffisant pour "n" tests Contains sufficient for "n" tests Contenuto sufficiente per "n" saggi Contém o suficiente para "n" testes		DE ES FR GB IT PT	Inhalt Contenido del estuche Contenu du coffret Contents of kit Contenuto del kit Conteúdo do kit
 Max Min	DE ES FR GB IT PT	Temperaturbereich Limitación de temperatura Limites de température de conservation Temperature limitation Limiti di temperatura Temperaturas limites de conservação		DE ES FR GB IT PT	Bestellnummer Número de catálogo Références du catalogue Catalogue number Numero di Catalogo Número do catálogo
	DE ES FR GB IT PT	Vor direkter sonneneinstrahlung schützen Mantener alejado de la luz solar Tenir à l'écart de la lumière du soleil Keep away from sunlight Tenere lontano dalla luce solare Mantenha longe da luz solar			

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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