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SURE Chemically Competent Cell 说明书

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E.coli B e14-(McrA-) Δ(mcrCB-hsdSMR-mrr)171 endA1 gyrA96 thi-1 supE44 relA1 lac recB recJ sbcC umuC::Tn5 (Kanr) uvrC [F' proAB lacIqZ Δ M15 Tn10 (Tetr)].

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Stratagene standard protocol

1. Pre-chill a 14-ml BD Falcon polypropylene round-bottom tubes on ice. Preheat NZY+ broth to 42°C .
2. Thaw the cells on ice. When thawed, gently mix and aliquot $100 \mu\text{l}$ of cells into the pre-chilled tubes.
3. Add $4 \mu\text{l}$ of the β -ME (β 巯基乙醇) to the aliquot of cells.
4. Swirl the tubes gently. Incubate the cells on ice for 10 minutes, swirling gently every 2 minutes.

5. Add 0.1-50 ng of the experimental DNA (or 2 μ l of a ligation mixture) to the aliquot of cells.
6. Swirl the tubes gently, then incubate the tubes on ice for 30 minutes.
7. Heat-pulse the tubes in a 42°C water bath for 30 seconds. The duration of the heat pulse is critical.
8. Incubate the tubes on ice for 2 minutes.
9. Add 0.9 ml of preheated (42°C) NZY+ broth and incubate the tubes at 37°C for 1 hour with shaking at 225-250 rpm.
10. Plate ≤ 200 μ l of the transformation mixture on LB agar plates containing the appropriate antibiotic (and containing IPTG and X-gal if color screening is desired).
11. Incubate the plates at 37°C overnight. If performing blue-white color screening, incubate the plates at 37°C for at least 17 hours to allow color development (color can be enhanced by subsequent incubation of the plates for 2 hours at 4°C).